


Biofinity[®] multifocal Fitting Guide



ADD	Dominant Eye	Non-Dominant Eye
+1.00D	D	D
+1.50D	D	D
+2.00D	D	N
+2.50D	D	N

Visual acuity expectations when using
D and N lens combination

Lens	Distance	Near
Binocularly	20/20	20/20
D Lens	20/20	20/40 or better
N Lens	20/40 or better	20/20

Clinical Tips

- Prescribe maximum plus power for distance vision (do not over minus)
- Choose the lower ADD power when possible; not necessary to overprescribe the ADD power
- Test patient's near function vision with their cell phone
- Check visual acuity with room lights on

3 Steps to Successfully Fit Biofinity® multifocal

- ① Start with a new refraction and verification of eye dominance (fogging technique).
- ② Select the distance prescription based on spherical equivalent corrected for the vertex distance. Choose D or N lens design based on needed ADD power (see table on reverse side).
- ③ Allow patients to adapt to lenses for a minimum of 15 minutes before assessing vision. If binocular vision is unacceptable, perform a monocular over refraction, using hand-held trial lenses, to determine which eye needs improvement.

To improve distance vision add $\pm 0.25\text{D}$ (up to $\pm 0.50\text{D}$) to the eye that needs improvement.

To improve near vision add $\pm 0.25\text{D}$ (up to $\pm 0.50\text{D}$) to the eye that needs improvement.



CooperVision®